



Instruction Manual

Model 106C/130/155 Optical Module

Berkeley Nucleonics Corporation
2955 Kerner Blvd. San Rafael, CA 94901
Ph: 415-453-9955 Fx: 415-453-9956
www.berkeleynucleonics.com

WARRANTY

Berkeley Nucleonics Corporation warrants all instruments, including component parts, to be free from defects in material and workmanship, under normal use and service for a period of one year. If repairs are required during the warranty period, contact the factory for component replacement or shipping instructions. Include serial number of the instrument. This warranty is void if the unit is repaired or altered by others than those authorized by Berkeley Nucleonics Corporation.

CONTENTS

	Page
SECTION 1 SPECIFICATIONS	7
<u>Model 106C, 155 and 130 Characteristics</u>	7
<u>Module Status Byte Summary</u>	9
SECTION 2 OPERATING INFORMATION	10
<u>Features</u>	10
<u>General</u>	11
Power Up	11
Module Installation	11
Safety Precautions	11
Warm Up Requirements	11
Troubleshooting	11
Default Settings	12
Front Panel Description	12
LED Indicators	12
Connectors	13
Rear Panel Description	13
Mainframe Operation	13
Front Panel Programming	13
Remote Programming	16
SECTION 3 THEORY OF OPERATION	17
<u>General</u>	17
Module Interface	17
Laser Driver Output	17
<u>Circuit Description</u>	18
Module Interface Board	18
Laser Output Board	20
Inputs and Control Signals	20
Laser Module	21
CW Mode	21
External Modulation Mode	23
Pulse Mode	24
Pulse Mode (Nonzero Baseline)	24
Pulse Mode (Zero Baseline)	25
External Drive	25
Impulse Mode	25
Laser Protection	26

CONTENTS

SECTION 4 MAINTENANCE AND CALIBRATION	27
<u>Maintenance</u>	27
Light Output Connector	27
<u>Calibration</u>	27
General	27
Equipment Required	27
Procedure	28
Module Interface DAC Calibration	28
External Drive Discriminator	28
CW and External Modulation Calibration	29
Pulse Baseline Calibration	29
Pulse Peak Calibration	29
Pulse Dynamic Characteristics	30
SECTION 5 PARTS LISTS AND SCHEMATICS	31
<u>Parts List</u>	31
Laser Output Board, 155-1	31
Module Interface Board, 6040-4	35
<u>Schematics</u>	<u>Dwg.No.</u>
Interconnection Diagram	155-30
Module Interface Board	6040-35A
Laser Output Board	155-31

CONTENTS

ILLUSTRATIONS

<u>Figure No.</u>		<u>Page</u>
1-1	Trigger and Output Pulse Timing	
2-1	Safety Labels	
3-1	Module Interface Block Diagram	
3-2	Laser Output Block Diagram	

TABLES

<u>Table No.</u>		<u>Page</u>
1-1	Module Status Byte Summary	9
2-1	Menu Keys for the 135/130 Module	14
3-1	Plug-In Module Memory Map	19
3-2	Control Signals	22

1064nm, 1560nm AND 1310nm OPTICAL MODULES

MODEL 106C, MODEL 155 AND MODEL 130

Graphic (Model 155 & 130)

The Model 106C, Model 155 and Model 130 are three in a series of plug-in modules that provide electrical and optical output pulses when installed in the Model 6040 mainframe.

These particular modules provide optical pulses of 1064, 1560 and 1310 nm wavelengths at peak levels to 1 mV at rates to 100 MHz.

SECTION 1

SPECIFICATIONS

MODEL 106C/155/130 CHARACTERISTICS

Timing Characteristics

Rep Rate:	0 Hz-100 MHz
Width:	3 ns - 640 s (Pulse Mode); 3 ns (min.) at reduced amplitude. Impulses, fixed 400 ps fwhm (typical).

Input Characteristics

EXTERNAL DRIVE

Range:	dc - 300 MHz (200 MHz for zero Baseline level): specifications apply dc - 100 MHz.
Input Impedance:	50 ohms.
Minimum Signal Amplitude:	300 mV.
Maximum Signal Amplitude:	± 7 V dc or 10 V ac p-p.
Minimum Width:	5 ns; 3 ns (min.) at reduced amplitude.
Threshold Range:	± 2.5 V.
Threshold Resolution:	10 mV.
Insertion Delay:	5 ns, typical (between EXT DR and LIGHT OUT).
Jitter:	30 ps rms (between EXT DR and LIGHT OUT).
Connector:	Threaded SMA (3 mm).

EXTERNAL MODE

Range:	100 Hz to 700 MHz (-3 dB).
Input Impedance:	50 ohms.
Sensitivity:	200 mV rms/mW, typical.
Maximum Signal Amplitude:	± 2 V dc or 3 V ac p-p.
Total Harmonic Distortion:	< 20 dB below fundamental (0.5 mW avg, level with 0.5 mW rms modulation).
Insertion Delay:	5 ns, typical (between EXT MOD and LIGHT OUT).
Connector:	Threaded SMA (3 mm).

SPECIFICATIONS cont'd.

Output Characteristics

LIGHT OUT

Wavelength:	Model 106C: 1064nm \pm 30 nm Model 155:1550 nm \pm 30 nm. Model 130:1310 nm \pm 50 nm.
Spectral Width:	2 nm rms (from 50 μ W to 1 mW).
Power Level:	1 mW max. (Peak or Baseline). 0 mW min. Impulse is fixed at 50 μ W Baseline and 0.5 mW Peak (typical).
Power Level Resolution:	5 μ W (from 50 μ W to 1 mW),
Extinction Ratio:	CW Mode and External Modulation Mode, has 5 μ W resolution from 0 to 1 mW
Mesial Level:	∞ (zero Baseline), or 20.0 to 1.01.
Pulse Adjustment Range:	25 μ W (zero Baseline), or 55 μ W to 0.995 mW. (Max. pulse size: min. pulse size.) 20 dB (non-zero Baseline), 13 dB (zero Baseline), 20 dB (CW. External Modulation).
Accuracy, Absolute:	\pm 1 dB (from 50 μ W to 1 mW).
Accuracy, Relative:	\pm 0.5 dB (\pm 10 %) from 100 μ W to 1 mW (Pulse Mode). \pm 0.05 dB (\pm 0.5 %) from 50 μ W to 1 mW (CW. External Modulation).
Temperature Coefficient:	0.05 dB/ $^{\circ}$ C.
Transition Times (10 to 90%):	0.5 ns rise time, 1 ns fall time (zero Baseline), 1 ns rise time, 1 ns fall time (nonzero Baseline).
Insertion Delay:	35 ns, typical (between mainframe TRIG OUT and LIGHT OUT; see Figure 1-1).
Jitter:	100 ps rms (between mainframe TRIG OUT and LIGHTOUT).
Connector:	ST, 8/125 μ m single-mode fiber, 0.12 NA (contact factory for other connectors).

SPECIFICATIONS cont'd.

Modes

PULSE

Conventional pulse generator with rate, delay, width and single/double pulse selections controlled by the 6040 mainframe. External Drive operation produces pulses corresponding in rate and duty cycle to an external pulse train.

EXTERNAL MODULATION

Converts digital and analog electrical signals are into their optical equivalent.

IMPULSE

Provides a narrow pulse of fixed width and amplitude, with rate and delay controlled by the 6040 mainframe.

Module Status Byte Summary

Table 1-1. Module Status Byte

<u>Bit</u>	<u>Description</u>
7	Always zero
6	Always zero
3	Always zero
4	Always zero
3	Always zero
2	Always zero
1	Laser Active
0	Laser Guard

SECTION 2

OPERATING INFORMATION

FEATURES

The Model 106C, 155 and 130 plug-in modules provide 1064nm, 1550 nm and 1310 nm optical output sources for the Model 6040 Universal Pulse Generator. Accurate and adjustable outputs are available for all of the four modes in which the 6040 mainframe can operate.

In Pulse Mode operation, the 106C/155/130 supplies flat-topped pulses with fast rise and fall times and independently adjustable Peak and Baseline levels. The timing for these pulses may be supplied in a number of ways.

A delayed pulse of adjustable width may be generated by the mainframe. This delay can be specified with respect to an internal trigger occurring at a selected repetition rate or with respect to an externally supplied trigger signal (TRIG IN). In addition, Single Cycle operation allows the user to trigger the instrument manually, using a pushbutton (or using remote programming). In each case a trigger out signal (TRIG OUT) is provided by the 6040 for reference (see Figure 1-1). Double Pulse operation, producing both an initial and delayed pulse out of the same jack may be selected for any trigger choice. For all of these timing options, the mainframe supplies electrical output pulses (PULSE OUT) coincident with the module's optical output (LIGHT OUT).

External Drive operation is also available in the Pulse Mode. This allows a drive signal, supplied by the user to the module's front panel (EXT DR), to generate the optical pulses directly. The occurrence and duration of each light pulse will correspond to that of each pulse in the external drive signal: when the EXT DR pulse goes high, the light pulse goes to the Peak value. Peak and Baseline output levels remain specified by the mainframe (as well as the threshold level for the external drive signal). With External Drive, the mainframe's pulse and trigger outputs are disabled.

In External Modulation Mode, an electrical signal supplied to the module's front panel (EXT MOD) by the user will be converted into its optical analog. The quiescent optical level, corresponding to the zero point that the input signal modulates about, is selected by the mainframe (External Modulation Level), but no other parameters can be altered.

In the Impulse Mode, a narrow optical pulse of fixed width and amplitude is produced at the module output (LIGHT OUT) with a corresponding pulse (of 5 ns duration) appearing at the mainframe (PULSE OUT). As with the Pulse Mode, the trigger source may be internal, external or from manual (or remote programming) control. Either a single delayed impulse or a pair of impulses separated by a delay may be obtained. In Impulse Mode, External Drive operation is disabled; and Peak and Baseline Level settings also have no effect.

CW Mode results in a steady-state optical output. The output power level in this Mode may be adjusted by the mainframe (CW Level).

For more detailed information on the characteristics of each Mode and how to control the module from the mainframe, see the Specifications section and the 6040 manual.

OPERATING INFORMATION

General

POWER UP

When power is applied to the 6040 mainframe with a 106C, 155 or 130 module installed, the instrument settings from the module's memory 0 are activated. The mainframe automatically checks what type of plug-in module is in place and loads the appropriate parameters. The LCD, after showing the mainframe's software version number and performing a memory check, will display "106C Ver. x.x" (or "155 Ver. X.x" or "130 Ver. x.x") where "x.x" is the version number of the module.

Module Installation

The module must be installed with mainframe power off. A module can be damaged or have its memory corrupted if inserted or removed from the mainframe with the power on. To install the module, simply slide it in and tighten the mount screw knob.

Safety Precautions

Laser light emitted from the end of the connected light fiber is invisible. Fibers should be terminated in a system which will not allow human exposure to this radiation. Do not stare into the beam or into a beam from a reflecting surface. Use of controls or adjustments or performance procedures other than those specified herein may result in hazardous exposure. Safety labels attached to the module are shown in Figure 2-1.

Always keep the LIGHT OUT connector covered with the dust cap when a fiber is not attached in order to avoid hazardous exposure and to keep the connector as clean as possible.

Warm Up Requirements

No warm up period, other than that necessary for the mainframe itself, is required.

TROUBLESHOOTING

Follow the procedure in the Troubleshooting section of the 6040 manual. Make sure that the module is seated correctly in the mainframe. If the module's memory has been corrupted, a cold boot will put the module's default settings (given below) into effect (as described in the Cold Boot paragraph in the 6040 manual).

OPERATING INFORMATION cont'd

The Quick Test procedure for the mainframe may be applied to the 106C/155/130 by selecting the Pulse Mode and following the test sequence using the module's LIGHT OUT connector and an optical detector in place of the mainframe's PULSE OUT. After Pulse Mode operation has been verified, the Impulse Mode can be tested. This will require a sampling oscilloscope with a 1 GHz bandwidth. The impulse should appear as a narrow pulse (approximately 400 ps). Its amplitude is roughly equal to that of a 0.3 mW, 5 ns pulse generated in the Pulse Mode. The CW and External Modulation Modes are similar. These should produce continuous outputs. The output under External Modulation should follow any waveform presented to the module's EXT MOD connector.

Default Settings

The following default settings for the 106C/153/130 module go into effect whenever a cold boot is performed.

MODE: Pulse

TRIG: Single Cycle
(with other values set as follows)
Internal Trigger Rate = 1 kHz.
External Trigger Threshold = 0.00 V
Trigger Slope +
External Drive Threshold = 0.00 V

TIMING: Delay = 1 μ s
Width = 1 μ s
Single Pulse

LEVEL: Peak= 1.000 mW
Baseline = 0.000 mW
External Modulation
CW=0.00 mW

GPIB/RS232: GPIB Address = 6
Baud Rate = 1200
Full Duplex
Remote Enabled

Front Panel Descriptions

LED INDICATORS

The 106C/155/130 module has two LED indicators on the front panel.

LASER ACTIVE will light when the instrument is in a state that can produce a light output.

LASER GUARD will light if the laser is being protected from an unsafe power output level, whenever the Mode is changed and when the unit is first turned on. LASER GUARD will remain on if any unsafe circuit conditions exist (e.g., excessive current driving the laser).

OPERATING INFORMATION cont'd

CONNECTORS

Three connectors appear on the 106C/155/130 module front panel.

LIGHT OUT provides the optical output from the module. A single-mode 8/123 Mm fiber with an ST connector is required (the unit can be configured for other connectors; consult the factory for details). To keep dust out of the connector when a fiber is not attached, a dust cap is provided.

EXT DR (External Drive) is a threaded SMA (3 mm) connector for accepting drive signal inputs to the module.

EXT MOD (External Modulation) is a threaded SMA (3 mm) connector that accepts external signals to modulate the instantaneous light output power level.

Rear Panel Description

The rear panel of the 106C/155/130 module has a mounting screw (for installation into the mainframe), one 40-pin edge connector, and one SMA connector.

The 40-pin connector allows the 6040 mainframe to control and communicate with the module and also supplies the power to the module.

The SMA (snap-on type) connector receives the high speed pulse generator DRIVE signal from the mainframe. This signal is an ECL version of the mainframe's front panel PULSE OUT.

Mainframe Operation

This section presents information on how to operate the Model 6040 Universal Pulse Generator with the Model 106C, 155 or Model 130 optical module installed. Only the details that are specific to the module are described. For an overall description of how to use the mainframe with plug-in modules, please refer to the 6040 manual.

Front Panel Programming

When operating the instrument from the front panel, certain control keys have module dependent action, as indicated in the 6040 manual. The aspects of these keys, and the menus they control, that are not general to the mainframe will be listed here.

An overall chart of the menu keys, showing which menu selections have control in each Mode, is given in Table 2-1. An "x" in the column for a given Mode indicates that the menu selection operates in that Mode.

OPERATING INFORMATION cont'd

Table 2-1. Menu Keys for the 133/130 Module.

<u>MODE Menu</u>				
	Pulse	Impulse	CW	External Modulation
TRIG Menu Single Cycle Internal Trigger (and Rate) External Trigger (and Threshold) External Trigger Slope External Drive (and Threshold)		X X X X		
TIMING Menu Delay Width Single/Double Pulse	X X X	X X		
LEVEL Menu Peak Baseline CW External Modulation	X X		X	X

OPERATING INFORMATION cont'd

{MODE} The Mode menu for the 106C/155/130 has all four selections available: Pulse, Impulse, CW and External Modulation.

Pulse Mode can operate over the entire timing range of the 6040, producing flat-topped delayed pulses. The Delay interval, Peak level, Baseline level and pulse Width are all adjustable.

External Drive operation is available in this Mode, allowing the module to be digitally modulated at rates from zero to 300 MHz (for nonzero Baseline levels) or 200 MHz (for zero Baseline level).

Impulse Mode produces pulses of fixed 400 ps width with fixed 30 uV baseline and 1 mW peak levels (typical).

CW Mode provides continuous wave optical output.

External Modulation Mode allows a user-provided analog or digital signal to be linearly converted to its optical equivalent. Inputs ranging from 100 Hz to 700 MHz are converted with 200 mV rms/mW sensitivity (typical).

{TRIG} The Trigger source and parameter menu, which operates in the Pulse and Impulse Modes, has all five menu items: Single/Double Cycle, Internal Trigger (and Rate), External Trigger (and Threshold), External Trigger Slope, and External Drive (and Threshold).

Internal Trigger Rate and External Trigger Threshold are adjustable over the 6040's entire range.

External Drive is valid only in Pulse Mode (and has no effect in Impulse Mode). The External Drive Threshold for the discriminator on the EXT DR input to the module is adjustable from -2.3 V to +2.3 V with 10 mV resolution.

{TIMING} This module places no constraints on timing settings selected with the Timing parameter menu.

OPERATING INFORMATION cont'd

{LEVEL}

All four Level parameter menu selections are available for this module: Peak Level, Baseline Level, External Modulation Level and CW Level.

Peak and Baseline Levels select the high and low power levels in Pulse Mode operation. These levels may be set to zero or they may be adjusted between 30 μW and 1 mW in 3 MW steps (levels between zero and 50 μW may be selected but output characteristics are not guaranteed). Peak Level and Baseline Level may be set independently, but if Peak is set below Baseline the output will be a constant CW at the Baseline Level.

External Modulation Level sets the quiescent optical power level. This may be selected over the entire range of zero to 1 mW in 3 MW steps.

CW Level may also be adjusted to any power level between zero and 1 mW with 5 μW resolution.

{UNITS}

This key is not used with the 106C/155/130 module (and has no effect).

FUNCTION KEYS

These keys are not used with the 106C/155/130 module (and have no effect).

{A}, {B}, {C}

Remote Programming

The remote programming commands that are specific to this module correspond to the module dependent front panel commands, as described in the previous section. Consult the 6040 manual and the above section for details on controlling the instrument with remote programming. The Module Status command, the only command that is specific to the 155/130 module, will be described here.

PS

Module Status: This command returns the Module Status byte for the 155/130.

Module Status Byte:

<u>Bit</u>	<u>Description</u>
7	Always zero
6	Always zero
5	Always zero
4	Always zero
3	Always zero
2	Always zero
1	Laser Active
0	Laser Guard

Bits 2-7: These bits are always zero and are reserved for future use.

Bit 1: This bit is set if the LASER ACTIVE LED is lit.

Bit 0: This bit is set if the LASER GUARD LED is lit.

SECTION 3

THEORY OF OPERATION

GENERAL

Module Interface

Figure 3-1 shows a simplified block diagram of the Module Interface board. The path for communication between mainframe and module is via PS. The eight QAD lines and five QA lines are the bus interface lines, and a MOD DIS line is used for disabling the Output board. Power is also delivered by P8.

The address demultiplexer and select logic circuitry decodes the bus signals and selects one of the other blocks.

The I.D. ROM contains information necessary for operation specific to the Output board. This includes boundaries for parameters, values used to initialize the nonvolatile RAM (NVRAM), and the version number of the I.D. ROM.

The nonvolatile RAM is used for saving and retrieving ten panel settings. It also holds the GPIB/RS232 bus settings.

The digital control circuits are used to monitor and set the operating state of the Output board. The DACs and amplifiers allow setting of up to four analog values used for level control. All of the control, status, and analog signals are delivered to the Output board by two 20-pin connectors, J1 and J2.

Laser Driver Output

Figure 3-2 shows a simplified block diagram of the Laser Output board. The connection to the Module Interface board is shown in the lower left hand corner. There are three analog lines used to set the output levels and the threshold of the External Drive discriminator. The control lines select which Mode the Laser Output will operate in.

The selection of the digital drive source is controlled by three digital control lines, EXT DRIVE ENABLE, IMPULSE, and EXT MOD. These determine which of the two drive sources will be controlling the state of S1, (a high speed transistor switch). A discriminator, whose threshold is set by the analog voltage EXT DRIVE LEVEL, is connected to the drive source from the front panel EXT DR connector. The other signal source is the rear panel DRIVE connector which delivers the mainframe's pulse generator output.

The level control circuitry sets one or both of the current sources as determined by the analog inputs PEAK LEVEL, and BASELINE LEVEL and the Mode control lines. During the Impulse, CY and External Modulation Modes, it uses the PIN detector in the laser module to monitor and stabilize the laser's output. When LASER GUARD is active, the current sources are set for zero output.

THEORY OF OPERATION

The two current sources are used for different Modes. The current source on the left is used during Pulse and Impulse Modes (to supply the Peak level), and during the External Modulation Mode. This current source can be modulated from a wideband preamp which is driven by the front panel SMA connector EXT MOD. The current is switched to the laser by SI, as determined by the DIGITAL DRIVE SELECT and the drive signals. The current source on the right is used during PULSE MODE when a nonzero Baseline level is selected and in the CW Mode.

The Laser Protection circuitry monitors the laser output and circuit conditions. If the laser output exceeds a preset level, S2 will be closed and all laser drive current will be shunted to ground. The current is then monitored until it is reduced below another preset value before S2 is allowed to reopen. This block drives the front panel LEDs, LASER ACTIVE and LASER GUARD, which allow the user to monitor the state of these circuits. The DAC EN signal is used to notify the Module Interface board of a potential problem and disable its DACs, thus setting all level control voltages to zero. The MOD DIS signal is used to disable the module's output.

The Laser Module is a single package that contains the laser diode, a PIN monitor diode, a thermoelectric cooler (TEC), and a temperature sensing device. An 8/123 Um single-mode fiber is connected to the front panel (LIGHT OUT). The TEC and thermistor are used in conjunction with an op amp and transistor to regulate the temperature of the laser diode, thus improving power and wavelength stability.

CIRCUIT DESCRIPTION

Module Interface Board (Schematic 6040-35)

The Module Interface board contains all the necessary circuits to allow the 6040 mainframe to control the module. Interfacing between the module and the mainframe is realized via the 40-pin edge connector, P5. This delivers eight data bits (QADO-QAD7) and 13 address bits (QADO-QAD7, multiplexed, and QA8-QA12, nonmultiplexed). P8 also delivers power to the Interface board and the Output board.

Z2 is an eight bit latch that demultiplexes QADO-QAD7 to produce the lower eight address bits QAO-QA7. This allows up to 8K bytes to be addressed within the Interface board, though not all of this is utilized at the present time. Table 3-1 gives the memory map for the module.

Z1 is a dual quad selector. Z1 A is used to select between one of four 2K byte segments. These four segments are used for: 1) Z3, the I.D. ROM; 2) Z4, the nonvolatile RAM; 3) reserved for expansion; 4) Z5, Z6 and Z7, the digital and analog control of the output board. Z1B selects which IC in the fourth segment is accessed.

Z3 is the I.D. ROM. It is an 8K byte ROM (only the lower 2K is used) that contains the module dependent information necessary for the 60-40 to operate correctly.

THEORY OF OPERATION

Table 3-1. Plug-In Module Memory Map

Memory Range

C000-C777	Z3, I.D.ROM
C800-CFFF	Unused
D000-D7FF	Z4, Nonvolatile RAM (NVRAM)
D800-DFFF	I/O
D800-DFFF	
D800-D9FF	Z5, 82C55 PPI
D800	Port A
D801	Port B
D802	Port C
D803	Control
DA00-DBFF	Unused
DC00-DDFF	Z6, 7528; Dual DAC
DC00-DC03	PEAK LEVEL
DC04-DC07	BASELINE LEVEL
DE00-DFFF	B27, 7528; Dual DAC
DE00-DE03	SPARE
DE04-DE07	EXT DRIVE LEVEL

THEORY OF OPERATION

This includes the Modes that are valid, parameter boundaries, and the type of output that the module has (optical or electrical). It also contains the values for initializing the nonvolatile RAM.

Z4 is a 2K byte nonvolatile RAM (NVRAM). It is used to save instrument settings and power-on conditions.

Z5 is a configurable Parallel Peripheral Interface 1C set up to allow 16 bits of output and eight bits of input. The outputs are used to select the operating state of the Output board while the inputs monitor the board's status.

DACs Z6 and Z7, in conjunction with op amps Z8 and Z9, generate analog signals that set the value of the Peak, Baseline, and quiescent level of the output. They are also used to set the threshold level for the External Drive discriminator. The variable resistors, R3-R6, are used to compensate for slight differences in the DACs.

The voltage reference for the analog signals is CR1, a 6.2 V temperature compensated zener diode. Q1 is used to disable CR1 if the signal DAC EN is allowed to go up to +12 V. The Output board asserts this in the event of a Mode change or if any unusual condition is detected.

The MOD DISABLE line allows the Output board to be disabled via the 6040 rear panel MODULE DISABLE connector.

The Module Interface board has two grounds, digital and analog. The analog ground is used exclusively by the DAC and op amp circuitry (Z6, Z7, Z8, Z9, etc.), while all other circuitry is connected to digital ground. These are separated to prevent any noise or dc offsets by power supply return currents from affecting the analog control voltages.

The 20-pin connectors J1 and J2 deliver the digital and analog control signals to the Output board, while the 16-pin connector J11 delivers the power supply voltages. J3 is an expansion connector to be used in conjunction with future modules.

LASER OUTPUT BOARD (Schematic 155-31)

Inputs and Control Signals (Schematic Sheets 1 and 2)

The Laser Output board receives three types of inputs: digital control signals from the Interface board via P1, analog control signals from the Interface board via P2, and drive and external modulation inputs via front or rear panel coax connectors. Four of the digital control signals, IMPULSE, EXT MOD, CW, and BIAS, are used to control analog switches, Z14 and Z11 (sheet 1). These switches are in their low impedance (closed) state when the control input, the terminal with inversion circle, is low (near ground).

THEORY OF OPERATION

Table 3-2 shows how these digital control signals affect the condition of the analog switches in each Mode. The switches are identified by their control terminals (e.g., Z14-9); "L" and "H" indicate low and high logic/voltage levels. By turning on and off these switches, these digital signals provide proper routing for BASELINE LEVEL and PEAK LEVEL, the analog signals that control the amplitude at the laser output. As an example, in the Pulse Mode with a nonzero Baseline (and not using External Drive), the only high input to Z9 is BIAS (Z9-5). Thus Z9-10 and Z9-12 are both high and Z9-6 is low. Z14-1 is low and this connects the BASELINE LEVEL amplitude control to Z8-3 via Z13-7 and Z14-3 Other Modes may be similarly analyzed using this table.

Three digital control signals, EXT MOD, CW (sheet 1) and EXT DRIVE EN (sheet 2), select which drive source is presented at the output of the ECL multiplexer, Z4. In Pulse Mode either the DRIVE signal from the mainframe or the EXT DR signal from the module front panel can be selected. This signal, after passing through Z3 and R162, becomes PREDRIVE, which determines when the laser is switched between Peak and Baseline levels. In External Modulation Mode, EXT MOD selects X4 (Z4-11), which is tied high, for the multiplexer output, causing PREDRIVE to be held high. In CW Mode, X2 is selected, causing PREDRIVE to be held low.

Laser Module

The Laser Module is shown enclosed by dashed lines within which are six components. The laser itself is the left hand diode whose anode is grounded and connected to the case at pins 13 and 13. The case is connected to ground via pins 1 and 2. A 170 μ H choke is provided internally for applying dc bias to the laser. This is used in all the Pulse Modes with nonzero Baseline and in the CW Mode. The diode shown between pins 7 and 8 is a built in monitor (detector) that is used for controlling the laser's output in the CW, External Modulation, and Impulse Modes and is also used to detect excessive optical power. An internal 22 ohm resistor provides wideband termination for the high speed signals from 06. A thermoelectric cooler (TEC) is incorporated in the module and is shown connected between pins 3/4 and 3/6. Its associated temperature sensor is connected between pins 11 and 12. The temperature is maintained at approximately 20° C by means of a feedback loop consisting of Z1 and 010. This loop supplies enough cooling current to bring the sensor voltage to 2.9 V (Z1-3). R104 limits the maximum current during start-up to a safe level (approximately 2 A).

CW Mode (Schematic Sheets 1 and 3)

CW operation utilizes the dc current source, Q1 and Z7. The voltage from Z8-1 (labeled CW OR BASELINE LEVEL) is applied to one arm of the resistor bridge, R36, R37, R63 and R66. This same voltage is forced to appear across R82 and R83. A typical calibration is for R82 and R83 to be 30 ohms so that 20 mA/V is produced by Q1. This current flows (via R81, L1, and the internal 170 μ H inductor) through the laser to ground. Since PREDRIVE is held low the path for current through Q6 will be cut off. Once the threshold current of the laser is exceeded, the relationship between the current in the laser and its optical output is linear.

THEORY OF OPERATION

The path of the CW OR BASELINE LEVEL signal may be followed on schematic sheet 1. First, we note that the PEAK LEVEL control voltage is inverted by Z13-1 and delivered to Z16-3 via R20. Second, we determine the status of the switches that affect the CW OR BASELINE LEVEL signal. Z11-3, for example, is connected to three such switches. From Table 3-1 it is seen that its control signal is IMPULSE from Z9-2, and that it is high in the CW Mode. Thus, Z9-3 is not conducting. Z14-16 is similarly found to be low so this switch is conducting. Also, Z14-1 is high which renders this switch open. Since Z14-9 is low, the EXT MOD OR CW FEEDBACK signal from Z16-7 is applied to Z16-2. The output from Z13-1 is also applied to Z16-3. The output from Z13-1 is not influenced by R22 because Z14-6 is open and no voltage is applied to R22.

Table 3-2. Control Signals

IC	Pulse		Mode			Destination	Signal Label
	Zero Bsln	Nonzero BSlN	CW	EXT MOD	Impulse		
Z9-1	L	L	L	L	H		IMPULSE
Z9-2	H	H	H	H	L	Z3 (Sheet 2)	$\overline{\text{IMPULSE}}$
Z9-11	L	L	H	L	H	Z4 (Sheet 2)	EXT MOD
Z9-10	H	H	H	L	H	Z10-5	$\overline{\text{EM}}$
Z9-13	L	L	H	L	L	Z4 (Sheet 2)	CW
Z9-12	H	H	L	H	H	Z10-4	$\overline{\text{CW}}$
Z9-5	L	H	L	L	L		BIAS
Z9-6	H	L	H	H	H	Z14-1	$\overline{\text{BIAS}}$
Z10-6	L	L	H	H	L	Z14-8	CW + EM
Z10-3	H	H	L	L	H	Z14-9	$\overline{\text{CW + EM}}$

THEORY OF OPERATION

Analog Switches

Z14-1	H	L	H	H	H	$\overline{\text{BIAS}}$
Z14-8	L	L	H	H	L	CW + EM
Z14-9	H	H	L	L	H	$\overline{\text{CW + EM}}$
Z14-16	H	H	L	H	H	$\overline{\text{CW}}$
Z11-1	H	H	H	H	L	$\overline{\text{IMPULSE}}$
Z11-8	H	H	H	H	L	$\overline{\text{IMPULSE}}$

In summary, the PEAK LEVEL control voltage is inverted about ground by Z13-1 and also undergoes a gain reduction of six (-6 V from PEAK LEVEL becomes +1 V at Z13-1). The current from the monitor diode is converted at Z16-7 to approximately 1 V for a 1 mW optical output and is applied to Z16-2 where it will almost cancel the signal from Z13-1. If this does not happen, a large error signal appears at Z16-1 where it passes via Z14-4 and Z8-1 to the current generator, Z7 and Q1 (sheet 3)

External Modulation Mode (Schematic Sheets 1 and 3)

We have traced in detail, for CW Mode, the manner in which the PEAK LEVEL control voltage causes a predictable current to be applied to the laser. Similar considerations show how the quiescent light level for the External Modulation Mode is derived. The PEAK LEVEL control voltage is still used to set the desired optical level (around which the modulation will occur). Feedback via Z16-7 is used as with CW. Z1 1-1 is high (as for CW) but now Z14-16 is also high. This has the effect of routing the error signal (Z16-1) through Z8-7 (a precision rectifier).

It is then applied (sheet 3) as the ETT MOD OR PULSE AMPLITUDE signal to Z1-5 (a x2 attenuator followed by an adjustable x3 gain stage). Since PREDRIVE is held high, the laser's dc bias is now provided by Q9 (via 06) R108 samples the current which is compared to the voltage from Z1-7 by a bridge arrangement (R72, R73, R74, R95). Z2-7 generates the error signal applied to Q8 and then to Q9.

The actual modulation is applied from J801, via FET current amplifier Q11, to Q8 where it generates a base voltage for Q9 which (via R112) generates a predictable signal current through R108 and Q6 that finally drives the laser. Although this circuit appears to be direct coupled from J801, it acts as if it were ac coupled due to the action of Z2, which monitors the total laser current (flowing through R108). The low frequency cut off is approximately 100 Hz.

THEORY OF OPERATION

Pulse Mode

In Pulse Mode, two conditions exist: the circuitry involved when the Baseline level is set to zero is different from the circuits used with a nonzero Baseline. In Pulse Mode with zero Baseline, there is no optical output between pulses (during a logical "zero" there is zero light output). The upper level is set by the PEAK LEVEL control voltage which, in turn, determines the amount of current switched into the laser. In Pulse Mode with a nonzero Baseline, both the upper ("one") and lower ("zero") logic levels must be adjustable according to the Peak and Baseline settings.

Pulse Mode (Nonzero Baseline)

Pulse Mode with nonzero Baseline requires 1) producing a current that determines the optical power radiated between pulses (the Baseline), and 2) calculating the difference between the desired Peak power and the Baseline power, and then generating a pulse of corresponding amplitude. This pulse is added to the already present Baseline level to produce the desired Peak level. The Baseline current is generated and applied to the laser in the same manner as for CW operation, except that optical feedback is not used.

The path from the BASELINE LEVEL control voltage is through Z13-7, Z14-3 (conducting), past Z14-14 and Z1 1-3 (both open), to Z8-1 (unity gain) and finally to Z7 (sheet 3). As stated above, Z7 and Q1 form a current source, and this is what supplies the Baseline current to the laser.

The path from the PEAK LEVEL control voltage is through Z13-1 (sheet 1), Z16-1, past Z14-15 (open), to Z5-7, through CR7, Z1-7 (sheet 3), and finally to Z2-7. As stated above, Z2-7, Q8 and Q9 form a current source which also drives the laser. Pulse Mode differs from External Modulation in that the current is switched by Q6 and Q7 (as defined by the level of the PREDRIVE signal) and there is no optical feedback.

The pulse current (through Q6 and Q7) is determined for PEAK LEVEL by the precision rectifier formed by Z5 and CR7 and the two resistor dividers R31/R40 and R32/R39 (sheet 1). This circuit produces a voltage which is the difference between the voltages at Z16-1 (the scaled version of the PEAK LEVEL control voltage) and Z8-1 (the scaled version of the BASELINE LEVEL control voltage), yet not less than zero. The output, at the cathode of CR7, is applied to Z1-5 (sheet 3) which provides an adjustable gain of from $\times 0.5$ to $\times 2$. At this point, a 3 V difference between the control voltages has become 500 mV at Z1-7. Z2-7 and R108 (47 ohm) convert this to approximately 10 mA, which is the value needed for a 0.30 mW optical power change. Thus a 3 V difference between the PEAK LEVEL and BASELINE LEVEL control voltages is converted to a 0.30 mW optical step. Since this step is added to the Baseline level already present (via Z7, Q1) the Peak and Baseline levels are determined by their respective control voltages, provided that the "gains" (voltage to current conversion factors) of the two channels are set to be equal.

THEORY OF OPERATION

PREDRIVE, the timing signal from the multiplexer (Z4), is applied (from Z3-2) to the base of predriver Q4. Q4 and Q3 are a switching pair whose current is controlled by Z2-1 and Q3. The predrive current through Q4 and Q3 increases with increasing optical output and, as the current increases, the main drivers Q6 and Q7 receive larger switching voltages. This is achieved by utilizing a current source comprised of Z2-1 and Q3, which tracks the EXT MOD OR PULSE AMPLITUDE control voltage at Z1-7. A small offset is introduced into this current source so that at very small pulse levels there is sufficient predrive current to switch Q6 and Q7.

Pulse Mode (Zero Baseline)

With Zero Baseline, only the switching transistors Q6 and Q7 (sheet 3) supply current to the laser. The PEAK LEVEL control voltage is processed directly through Z13-1, Z16-1 and Z8-7. Switches Z11-8, Z14-1, Z14-9 and Z14-6 are open and Z14-8 is closed. Thus Z13-1 is a x6 attenuator while Z16-1 is a x1 voltage follower. Since the switch Z14-1 is open, the BASELINE LEVEL control voltage is not presented to Z8-3, which will now be held at ground by R24. This allows Z8-7 (through CR7) to be a x1 follower of the voltage at Z16-1

This has two implications. 1) the CW OR BASELINE LEVEL current source is turned off, which prevents any quiescent light and 2) the PEAK LEVEL control voltage adjusts the pulsed current to values that produce a power level between zero and 1.0 mW (the same range that the CW OR BASELINE LEVEL current source uses with a nonzero Baseline Pulse Mode, as determined by the BASELINE LEVEL control voltage.

External Drive (Schematic Sheet 2)

External Drive is used in the Pulse Mode to shift the source of the drive signal from the mainframe to the front panel EXT DR connector. The circuits behave the same as under ordinary Pulse Mode operation except that the PREDRIVE signal is now obtained (via the multiplexer Z4-5) from the External Drive circuit. Two control signals from P2, EXT DRIVE LEVEL and EXT DRIVE POL, select the trigger threshold and polarity for the signal presented to the EXT DR jack.

Impulse Mode

Impulse Mode is used to generate a narrow, fixed width and fixed amplitude pulse. The digital timing is obtained from PREDRIVE (sheet 2). The leading edge of the impulse signal Z3-2 is caused by DRIVE itself. R157 and Z3-15 provide a short time delay (less than 2 ns) that produces the trailing edge. The result is a short, positive pulse delivered to the base of Q4 (sheet 3).

The predrivers (Q4 and Q3) and the main drivers (Q6 and Q7) have both been programmed by PEAK LEVEL for a high (120% of full scale) current that represents the amplitude of the impulse waveform. In order to assure fast response, the laser is biased to just above its threshold current by an auxiliary feedback loop that holds the laser at this point regardless of duty factor variations. The optical output is monitored by internal detector diode DSI-7, and in Impulse Mode represents the average optical power. Both the Baseline level and the average power of the impulses themselves contribute to this signal. The latter contribution varies with rep rate (duty factor).

THEORY OF OPERATION

Since it is desired to stabilize the Baseline level only, a signal proportional to the duty factor is required. This signal is obtained from Z3-2 (via R163). Both signals are sent to the auxiliary feedback loop. The IMPULSE COMPENSATION and IMPULSE FEEDBACK signals (sheet 1) are combined in Z6 along with a dc level from R37 and are sent back (via Z11-3 and Z8-1) to the Baseline current source (Q1, sheet 2). R37 adjusts for the small (approximately 10 μ W) amount of residual light that is emitted between impulses. R42 adjusts the balance between the optical sample (from D51-7) and the digital sample (Z3-2) so that there is no change with varying rep rates.

Laser Protection (Schematic Sheets 1 and 3)

The laser is protected from excessive drive currents by a group of circuits referred to as LASER GUARD. These circuits ensure that whenever an unusual situation occurs the laser is promptly shunted by a low impedance. This action is called "crowbar." The path around the laser is from DS1-10 through LI, Q2, and R79 to ground. Even if both sets of driver transistors (Q6, Q9 and Q1) were to short circuit, the resulting current would be shunted around the laser by Q2.

If there is excessive current (20 mA or more) flowing in the crowbar transistor, Q2, it is sensed by Z7. This causes a positive voltage (via CR14) to reinforce the triggering of Z12. The circuit remains in the crowbar state until the over current condition is removed. Although the control circuits generally turn on "gracefully" without excessive current transients, a time constant is provided (C59, R101) that turns Q2 on when the power is initially applied.

Another input to Z12-4 is from the internal monitor diode DSI-7. This produces a voltage across R102 and R103 that is set to trigger Z12 when the optical output becomes too high.

When a change of Mode is made, the ENABLE signal (Z9-9, sheet 1) triggers the crowbar action. This is applied to a timer circuit, Z12, that acts to turn on the crowbar transistor Q2 (via 012) It also is triggered at initial power-on by CR-23 and R130 A CROWBAR signal is sent to Z17-2 (sheet 1). This is another timer that acts to further hold off the control circuitry and permit a smooth return to normal operation.

While Z17 is in its timing cycle, the LASER ACTIVE LED is turned off and a positive level from Z17-3 via Z15-1 turns on the LASER GUARD indicator. After several seconds the Z17 timer will recover (if no over current condition is present) and allow the control circuits to recover.

SECTION 4

MAINTENANCE AND CALIBRATION

MAINTAINENCE

Light Output Connector

For satisfactory performance, proper care in the use of the optical components is necessary. There must be no contamination to interfere with the passage of light through the fiber end connections and the bulkhead connector. The following procedures should be observed.

1. Minimize the number of times connecting and disconnecting the optical cable.
2. Keep the connectors absolutely clean.
3. Keep the output bulkhead connector and the user's connector capped when not in use.

CALIBRATION

General

The calibration of the 106C/155/130 module is in two parts: the first is for setting control voltages on the Module Interface board (PCB 6040-4); the second is for setting the Laser Output board (PCB 155-1) for the correct power output and impulse characteristics and to verify the External Modulation bandwidth. It is recommended that the calibration of the 155/130 module be verified every 12 months. The instrument requires no warm up period exceeding that of the 6040, for which ten minutes is suggested.

Equipment Required

- 3-1/2 digit (or better) DVM.
- BNC 6100 Optical Power Meter (or equivalent) calibrated for the appropriate wavelength.
- 1 GHz bandwidth oscilloscope (sampling or real time). A Tektronix 7000 Series with appropriate sampling plug-ins is satisfactory.
- 1 GHz bandwidth 1330 nm detector (InGaAsP). A Tektronix S-42 Optical Sampling Head can be used with the 7000 Series Oscilloscope.
- 1 GHz variable sine wave generator. An HP 8657A Signal Generator with a Type N to SMA adapter can be used.
- Threaded SMA (3 mm) terminated 30 ohm coaxial cables, 1 meter length.

MAINTENANCE AND CALIBRATION

- BNC terminated 30 ohm coaxial cables, 1 meter length.
- Variable dc voltage source (capable of ± 3 V into 30 ohms).
- Single-mode optical fiber patch cord terminated with appropriate connectors..

PROCEDURE

Note: This calibration should be carried out in the order presented.

Before starting, verify that the power supply voltages are at their nominal levels (+12 V ± 0.1 V, -12 V ± 0.1 V, +5 V ± 0.05 V, and -5-2 V ± 0.05 V) with the Module plugged in. Operate the Module on extension cables.

Module Interface DAC Calibration

Refer to Schematic 6040-35 (PCB 6040-4). Set the Mode to Pulse; set the Peak and Baseline levels both to 1.000 mW; set the Trigger to External Drive with a threshold of +2.50 V.

Connect the DVM between TPG (A GND) and Z9-1. Adjust R3 for a voltage of -5.000 V ± 2 mV.

Connect the DVM between TPG (A GND) and Z8-7. Adjust R4 for a voltage of -4.000 V ± 2 mV.

Connect the DVM between TPG (A GND) and Z8-7, Adjust R5 for a voltage of -4.000 V ± 2 mV.

External Drive Discriminator

All further measurements and adjustments refer to the Laser Output board. Schematic 155-31 (PCB 155-1).

Set the Mode to Pulse. Set Trigger to External Drive with a threshold of +0.99 V. Connect the DVM to Z18-1 (Schematic 155-31, sheet 2) and record the voltage. Now set the threshold for -0.99 V. Adjust R141 so that the DVM reading exactly complements the recorded value.

Connect a +1.00 V dc supply to the module EXT DR connector. Connect the DVM to Z5-8. Scan the TRIG External Drive threshold parameter from 0.90 V to +1.10 V and verify that Z5-8 changes from low to high (-1.7 V to -0.8 V) when the threshold voltage is set for 1.00 V ± 30 mV.

MAINTENANCE AND CALIBRATION

CW and External Modulation Calibration

Using the patch cord, connect LIGHT OUT to the 6100 Optical Power Meter. Set the 6100 for Average power measurement and the 0 dBm range.

Set the 6040 Mode to CW and the CW level to 1.000 mW. If necessary, adjust R11 (Schematic 133-31, sheet 1) to obtain a reading of 1.000 mW \pm 10 μ W on the 6100. Set the CW level for 100 μ W and verify the output to be within 10 μ W.

Set the 6040 Mode to External Modulation (with no input connected to EXT MOD) and set the External Modulation level for 1.000 mW. Verify an output of 1.000 mW \pm 10 μ W. Set the External Modulation level for 100 μ W and verify the output to be within 10 μ W.

Pulse Baseline Calibration

Set the 6040 Mode to Pulse, and the Trigger to Single Cycle.

- A) Set the 6040 Baseline level to 1.000 mW. Adjust R83 (Schematic 133-31, sheet 3) to obtain a reading of 1.000 mW \pm 20 μ W on the 6100.
- B) Set the Baseline level to 100 μ W. Adjust R22 as necessary to obtain 0.100 mW \pm 20 μ W.

Repeat steps A) and B) until both readings are within 20 μ W of their stated values.

Pulse Peak Calibration

- A) Set the Baseline level to 0.000 mW, Set the Trigger to External Drive with a threshold of -50 mV. Adjust R58 (Schematic 155-31, sheet 3) for an output of 1.000 mW \pm 20 μ W.
- B) Set the Baseline level to 100 μ W. Adjust R19 to obtain an output of 1.000 mW \pm 20 μ W.

Repeat steps A) and B) until both readings are within 20 μ W of their stated values.

Note: The following sections require high bandwidth paths; any improperly cabled terminated paths will yield incorrect results.

MAINTENANCE AND CALIBRATION

Pulse Dynamic Characteristics

Connect the 6040 TRIG OUT to the External Trigger input of the oscilloscope (bandwidth must be at least 1 GHz). Connect the 6040 LIGHT OUT to the detector (also 1 GHz). Connect the detector's output to channel A of the oscilloscope.

With the 6040 still in Pulse Mode, set the Peak level to 1.000 mW and the Baseline level to 100 μ W. Set TRIG for Internal Trigger, 100 kHz. Set Timing for a Width of 40 ns. and a Delay appropriate to view the detector output on the scope. Verify that leading and trailing edges have transition times (10-90%) \leq 1 ns.

SECTION 5

PARTS LIST AND SCHEMATICS

Abbreviations

CER	Ceramic	PF	Pico farad
COMP	Composition	SIP	Single Inline Package
DIP	Dual Inline Package	TAN	Tantalum
ELEC	Electrolytic	UH	Microhenry
FAC SEL	Value Set at Factory	UF	Microfarad
K	Kilohm	V	Working Volts
M	Megohm	VAR	Variable
MF	Metal Film	W	Watts
MIC	Mica	WW	Wire wound
MONO	Monolithic Ceramic		

----- NOTE -----

*The number in the **second column** is the
BERKELEY NUCLEONICS re-order number*

LASER OUTPUT BD. 155-1, MODEL 155 649-144
LASER OUTPUT BD. 155-1, MODEL 130 649-145

C1	122-002 1 μ F 10% 35 V TAN	C23	110-033 0.1μF \pm 20% 50V CER MONO
C2	110-019 0.5 μ F 20% 25 VCER	C24	122-015 33 μF 10% 35 V TAN
C3	122-009 0.68 μ F 10% 35 VTAN	C25	110-033 0.1 μF \pm 20% 50V CER MONO
C4	110-021 0.01 μ F 20% 18 V CER	C26	110-033 0.1 μF \pm 20% 50V CER MONO
C5	122-001 047 μ F 10% 35 V TAN	C27	110-033 0.1 μF \pm 20% 50V CER MONO
C6	110-021 0.01 μ F 20% 18 V CER	C28	110-033 0.1 μF \pm 20% 50V CER MONO
C7	NOT USED	C29	110-033 0.1 μF\pm 20% 50 V CER MONO
C8	NOT USED	C30	110-033 0.1 μF\pm 20% 50 V CER MONO
C9	NOT USED	C31	110-033 0.1 μF\pm 20% 50 V CER MONO
C10	NOT USED	C32	110-033 0.1 μF\pm 20% 50 V CER MONO
C11	120-015 33 μ F 10% 35 VTAN	C33	110-033 0.1 μF\pm 20% 50 V CER MONO
C12	110-033 0.1 μ F \pm 20% 50 V CER MONO	C34	110-033 0.1 μF\pm20% 50V CER MONO
C13	110-033 0.1 μ F \pm 20% 50 V CER MONO	C35	110-033 0.1 μF\pm20% 50 V CER MONO
C14	110-033 0.1 μ F \pm 20% 50 V CER MONO	C36	122-015 33 μF 10% 35 V TAN
C37	110-033 0.1 μ F \pm 20% 50 V CER MONO	C38	110-033 0.1 μF\pm20% 50 V CER MONO
C15	110-033 0.1 μ F \pm 20% 50 V CER MONO	C39	110-033 0.1 μF\pm20% 50 V CER MONO
C16	110-033 0.1 μ F \pm 20% 50 V CER MONO	C40	122-015 33μF 10% 35 VTAN
C17	110-033 0.1 μ F \pm 20% 50 V CER MONO	C41	110-033 0.1 μF\pm 20% 50 V CER MONO
C18	1100330.1 μ F \pm 20% 50 V CER MONO	C42	110-033 0.1 μF\pm 20% 50 V CER MONO
C19	110-033 0.1 μ F \pm 20% 50 V CER MONO	C43	110-033 0.1 μF\pm 20% 50 V CER MONO
C20	110-033 0.1 μ F \pm 20% 50 V CER MONO	C44	110-0330.1 μF\pm 20% 50 V CER MONO
C21	110-033 0.1 μ F \pm 20% 50 V CER MONO	C45	110-019 0.05 μF 20% 25 VCER
C22	110-033 0.1 μ F \pm 20% 50 V CER MONO		

PARTS LIST AND SCHEMATICS

C46	NOT USED	CR7	411-004 IN4152
C47	NOT USED	CR8	NOT USED
C48	NOT USED	CR9	NOT USED
C49	NOT USED	C10	412-009 IN52318
		C11	411-003 IN4005
		C12	411-004 IN4152
C50	NOT USED		
C51	NOT USED	CR13	411-004 IN4152
C52	110-019 0.05 μ F 20% 25 VCER	CR14	411-004 IN4152
C53	110-021 0.01 μ F 20% 16 VCER	CR15	411-004 IN4152
C54	110-033 0.01 μ F 20% 50 V CER MONO	CR16	411-004 IN4152
		CR17	412-009 IN52318
C55	110-033 0.01 μ F 20% 50 V CER MONO		
C56	110-011 0.001 μ F 10% 1KV CER	CR18	413-002 IN4738
C57	110-034 100PF \pm 10% 25 V CER	CR19	411-002 IN270
C58	110-019 0.05 μ F 20% 25 V CER	CR20	411-004 IN4152
		CR21	411-004 IN4152
C59	122-015 33 μ F 10% 35 V TAN	CR22	NOT USED
C60	110-034 100 PF \pm 10% 25 V CER		
C61	110-011 0.001 μ F 10% 1 KV CER	CR23	411-004 IN4152
C62	110-021 0.01 μ F 20% 16 V CER	CR24	411-004 IN4152
C63	110-019 0.05 μ F 20% 25 V CER	CR25	411-004 IN4152
C64	110-011 0.001 μ F 10% 1KV CER	DSI	665-011 MODEL 106LASER(1064NM)
C65	110-034 100 PF \pm 10% 25 V CER	DSI	665-010 MODEL 155LASER(1550NM)
C66	110-011 0.001 μ F 10% 1KV CER	DSI	665-009 MODEL 130LASER(1310 NM)
C67	110-019 0.05 μ F 20% 25 V CER		
C68	122-013 3.3 μ F 10% 15 V TAN	01	430-027 MP83848
		02	430-008 2N2905
C69	NOT USED	03	430-026 MP83640
C70	NOT USED	04	430-056 MM4409
C71	NOT USED	05	430-056 MM4409
C72	110-021 0.01 μ F 20% 16 V CER		
C73	110-034 100 PF \pm 10% 25 V CER	06	430-059 A500
		07	430-059 A500
C74	110-034 100 PF \pm 10% 25 V CER	08	430-056 MM4049
C75	NOT USED	09	430-059 A500
C76	NOT USED	010	430-060 2N6041
C77	NOT USED		
C78	NOT USED	011	431-006 SD210
		012	430-027 MPS3646
C79	NOT USED	013	430-027 MPS3646
C80	112-016 10PF 5% 500 V MICA	014	431-006 SD210
C81	112-003 47 PF 5% 500 V MICA	015	430-026 MPS3640
C82	112-021 5 PF 5% 500 V MICA		
C83	112-019 15 PF 5% 500 V MICA		
		R1	223-010 1KX 5 SIP RES NETWORK
CR1	NOT USED	R2	NOT USED
CR2	NOT USED	R3	NOT USED
CR3	412-009 IN52318	R4	NOT USED
CR4	411-004 IN4152	R5	NOT USED
CR5	412-012 IN5227		
		R6	222-059 29.4 K 1% 1/4 W MF
CR6	411-004 IN4152	R7	222-059 29.4 K 1% 1/4 W MF
		R8	NOT USED

PARTS LIST AND SCHEMATICS

R9	213-301 300 OHMS 5% 1/4 W COMP	R54	NOT USED
R10	213-332 33 K 5% 1/4 W COMP	R55	NOT USED
R11	244-035 2 K MULTITURN	R56	222-026 110 K 1% 1/4 W MF
R12	NOT USED	R57	222-026 110 K 1% 1/4 W MF
R13	NOT USED	R58	244-011 1 K MULTITURN
R14	NOT USED	R59	NOT USED
R15	NOT USED	R60	NOT USED
R16	NOT USED	R61	NOT USED
R17	222-022 4.99 K 1% 1/4 W MF	R62	NOT USED
R18	222-045 4.22 K 1% 1/4 W MF	R63	NOT USED
R19	244-035 2 K MULTITURN	R64	NOT USED
R20	213-103 10 K 5% 1/4 COMP	R65	222-026 110 K 1% 1/4 W MF
R21	213-473 47 K 5% 1/4 COMP	R66	222-026 110 K 1% 1/4 W MF
R22	244-036 10 K MULTITURN	R67	213-102 1 K 5% 1/4 W COMP
R23	NOT USED	R68	213-102 1 K 5% 1/4 W COMP
R24	213-104 100 K 5% 1/4 W COMP	R69	222-026 110 K 1% 1/4 W MF
R25	213-105 1 M 5% 1/4 W COMP	R70	222-055 102 K 1% 1/4 W MF
R26	213-102 1 K 5% 1/4 W COMP	R71	222-070 95.3 K 1% 1/4 W MF
R27	213-102 1 K 5% 1/4 W COMP	R72	222-026 110 K 1% 1/4 W MF
R28	213-104 100 K 5% 1/4 W COMP	R73	222-025 110 K 1% 1/4 W MF
R29	213-274 270 K 5% 1/4 W COMP	R74	222-026 110 K 1% 1/4 W MF
R30	213-131 130 OHM 5% 1/4 W COMP	R75	214-080 65 OHM 5% 1/8 W COMP
R31	222-071 1 M 1% 1/4 W MF	R76	222-051 10K 1% 1/4 W MF
R32	222-071 1 M 1% 1/4 W MF	R77	222-043 237 K 1% 1/4 W MF
R33	213-102 1 K 5% 1/4 W COMP	R78	222-051 10 K 1% 1/4 W MF
R34	213-105 1 M 5% 1/4 W COMP	R79	213-150 15 OHM 5% 1/4 W COMP
R35	213-105 1 M 5% 1/4 W COMP	R80	222-017 215 K 1% 1/4 W MF
R36	213-104 100 K 5% 1/4 W COMP	R81	213-820 82 OHM 5% 1/4 W COMP
R37	244-036 10 K MULTITURN	R82	213-330 33 OHM 5% 1/4 W COMP
R38	213-820 82 OHM 5% 1/4 W COMP	R83	244-032 50 OHM 20-TURN
R39	222-071 1 M 1% 1/4 W MF	R84	213-101 100 OHM 5% 1/4 W COMP
R40	222-071 1 M 1% 1/4 W MF	R85	222-026 110 K 1% 1/4 W MF
R41	213-513 51 K 5% 1/4 W COMP	R86	213-102 1 K 5% 1/4 W COMP
R42	244-036 10 K MULTITURN	R87	213-301 300 OHM 5% 1/4 W COMP
R43	213-512 5.1 K 5% 1/4 W COMP	R88	213-102 1 K 5% 1/4 W COMP
R44	213-201 200 OHM 5% 1/4 W COMP	R89	214-501 500 OHM 5% 1/8 W COMP
R45	213-391 390 OHM 5% 1/4 W COMP	R90	214-580 68 OHM 5% 1/8 W COMP
R46	213-512 501 K 5% 1/4 W COMP	R91	214-510 51 OHM 5% 1/8 W COMP
R47	213-102 1 K 5% 1/4 W COMP	R92	214-131 130 OHM 5% 1/8 W COMP
R48	213-472 4.7 K 5% 1/4 W COMP	R93	214-391 390 OHM 5% 1/8 W COMP
R49	NOT USED	R94	214-510 51 OHM 5% 1/8 W COMP
R50	222-014 499 OHM 1% 1/4 W MF	R95	222-026 110 K 1% 1/4 W MF
R51	222-014 499 OHM 1% 1/4 W MF	R96	213-106 10 M 5% 1/4 W COMP
R52	213-471 470 OHM 5% 1/4 W COMP	R97	213-102 1 K 5% 1/4 W COMP
R53	NOT USED	R98	214-510 51 OHM 5% 1/8 W COMP
		R99	213-102 1 K 5% 1/4 W COMP
		R100	213-103 10 K 5% 1/4 W COMP

PARTS LIST AND SCHEMATICS

R101	213-102 1K 5X 1/4 W COMP	R142	222-018 249K 1% 1/4 W MF
R102	213-102 1 K 5% 1/4 W COMP	R143	222-018 249 K 1% 1/4 W MF
R103	244-038 5 K MULTI TURN	R144	222-039 1 K 1% 1 / 4 W MF
R104	231-008 2 OHM 1% 10 WWW	R145	222-091 182 K 1% 1/4 W MF
R105	222-050 8.66 K 1% 1 /4 W MF		
		R146	213-103 10 K 5% 1 /4 W COMP
R106	213-103 10 K 5% 1 /4 W COMP	R147	214-391 390 OHM 5% 1/8 W COMP
R107	214-510 51 OHM 5% 1 /8 W COMP	R148	214-391 390 OHM 5% 1 /8/ W COMP
R108	214-470 47 OHM 5% 1 /8 W COMP	R149	213-510 51 OHM 5% 1 /4 W COMP
R109	214-470 47 OHM 5% 1/8 W COMP	R150	213-510 51 OHM 5% 1 /4 W COMP
R110	213-103 10 K 5% 1 /4 W COMP	R151	213-101 100 OHM5% 1 /4 W COMP
		R152	213-510 51 OHM 5% 1 /4 W COMP
R111	214-820 82 OHM 5% 1/8 W COMP	R153	213-102 1 K 5% 1 /4 W COMP
R112	214-220 22 OHM 5% 1/8 W COMP	R154	213-102 1 K 5% 1 /4 W COMP
R 113	214-100 10 OHM 5% 1/8 W COMP	R155	213-102 1 K 5% 1 /4 W COMP
R114	213-105 1 M 5% 1/4 W COMP	R158	213-102 1 K 5% 1 /4 W COMP
R 115	213-103 10 K 5% 1/4 W COMP		
R116	222-013 422 OHM 1% 1 /4 WMF	R157	213-101 100 CHM 5% 1 /4 W COMP
R117	213-511 510 OHM 5% 1 /4 W COMP	R158	214-131 130 OHM 5% 1/8 W COMP
R118	213-105 1 M 5% 1 /4 W COMP	R159	214-820 82 OHM 5% 1 /8 W COMP
R119	214-101 100 OHM 5% 1 /8 W COMP	R160	214-820 82 OHM 15% 1 /8 W COMP
R120	214-392 3.9 K 5% 1/8 W COMP	R161	214-131 130 OHM 5% 1 /8 W COMP
		R162	214-580 56 OHM 5% 1 /8 W COMP
R121	214-390 39 OHM 5% 1 /8 W COMP	R163	214-102 1 K 5% 1 /8 W COMP
R122	213-102 1 K 5% 1 /4 W COMP	R164	222-080 332 OHM 1% 1 /4 W MF
R123	213-222 22 K 5% 1 /4 W COMP	R165	214-820 82 OHM 5% 1/8 W COMP
R124	213-102 1 K 5% 1/ 4 W COMP	R166	214-131 130 OHM 5% 1 /8 W COMP
R125	222-042 2 K 1% 1/ 4 W MF		
		R167	214-391 390 OHM 5% 1 /8 W COMP
R126	222-019 274 K 1% 1/ 4 W MF	R168	214-391 390 OHM 5% 1/8 W COMP
R127	213-510 51 OHM 5% 1 /4 W COMP		
R128	213-105 1 M 5% 1/ 4 W COMP	21	440-168 LF412
R129	213-223 22 K 5% 1 /4 W COMP	22	440-168 LF412
R130	213-223 22 K 5% 1 /4 W COMP	23	440-134 1OHD102
		24	440-200 MC1OH 154P
		25	440-148 SPO685
R131	213-102 1 K 5% 1 /4 W COMP	26	440-168 LF412
R131	214-510 51 OHM 5% 1 /8 W COMP	27	440-168 LF412
R133	NOT USED	28	440-168 LF412
R134	NOT USED	29	440-161S74HCO4
R135	NOT USED	210	440-150SHC00
R136	NOT USED	211	440-170 DG201
R137	NOT USED	212	440-143 NE528
R138	222-051 10 K 1% 1 /4 W MF	213	440-168 LF412
R139	222-051 10 K 1% 1 /4 W MF	214	440-170 DG201
R140	222-018 249 K 1% 1 /4 W MF	215	440-168 LF412
		216	440-168 LF412
R141	244-034 200 OHM MULTI TURN	217	440-144 555
		218	440-165 LF412

PARTS LIST AND SCHEMATICS

MODULE INTERACE BD 6040-4

649-143

C1 122-016 10UF±10% 15 V TAN
 C2 122-016 10UF±10% 15 V TAN
 C3 122-014 33UF±10% 6 V TAN
 C4 122-014 33UF±10% 6 V TAN
 C5 110-033 0.1UF±20% 50 V CER MGNO
 C6 110-033 0.1 UF±20%50V CER MONO
 C7 110-033 0.1 UF±20%50V CER MONO

C8 110-033 0.1 UF±20%50V CER MONO
 C9 110-033 0.1 UF±20%50 V CER MONO
 C10 110-033 0.1 UF±20%50V CER MONO
 C11 110-033 0.1 UF±20%50V CER MONO
 C12 110-033 0.1 UF±20%50V CER MONO
 C13 110-033 0.1 UF±20%50V CER MONO

C14 110-033 0.1 UF± 20% 50 V CER MONO
 C15 110-033 0.1 UF± 20% 50 V CER MONO
 C16 122-014 33 UF ± 10% 6 V TAN
 C17 122-014 33 UF ± 10% 6 V TAN
 C18 110-011 0.001 UF±10% 1 KV CER

C19 110-011 0.001 UF± 10% 1 KV CER
 C20 110-011 0.001 UF± 10% 1 KV CER
 C21 110-011 0.1 UF± 20% 50 V CER

 CR1 412-001 1N821

 Q1 430-026 MPS3540

R1 223-010 4.7 KX9 SP RES NETWORK
 R2 223-016 4.7 KX9 SP RES NETWORK
 R3 244-011 1 K PC MT MULTITURN
 R4 244-011 1 KPC MT MULTITURN
 R5 244-011 1 KPC MT MULTITURN
 R6 244-011 1 KPC MT MULTITURN
 R7 222-018 240 K 1% 1 /4 W MF

R8 222-018 249 K 1% 1 /4 W MF
 R9 222-018 249 K 1% 1 /4 W MF
 R10 222-018 249 K 1% 1 /4 W MF

R11 NOT USED
 R12 NOT USED
 R13 NOT USED
 R14 NOT USED
 R15 213-431 430 OHM 5% 1 /4 W COMP
 R16 213-103 10 K 5% 1/4 W COMP

 21 440-212 74HC139
 22 440-175 74HC373
 23 440-210 2764
 24 440-190 MK48202-20
 25 440-195 P82C55-2

26 440-213 AD7528
 27 440-213 AD7528
 28 440-168 LF412
 29 440-168 LF412